



WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 2

Catalog No: #AP88833

Package Size: #AP88833-1 50ug #AP88833-2 100ug #AP88833-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Accession No.	Uniprot ID:Q08E66
Uniprot	Q08E66
GeneID	531979;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:WFIKKN2

Note: This product is for in vitro research use only