

Transmembrane protein 132D

Catalog No: #AP88837

Package Size: #AP88837-1 50ug #AP88837-2 100ug #AP88837-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transmembrane protein 132D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Mature oligodendrocytes transmembrane protein,Mature OL transmembrane protein,HBE120,KIAA1944,MOLT
Accession No.	Uniprot ID:Q14C87
Uniprot	Q14C87
GeneID	121256;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TMEM132D

Note: This product is for in vitro research use only