

Mitochondrial uncoupling protein 3

Catalog No: #AP88838

Package Size: #AP88838-1 50ug #AP88838-2 100ug #AP88838-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitochondrial uncoupling protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	UCP 3,Solute carrier family 25 member 9,Slc25a9
Accession No.	Uniprot ID:P56501
Uniprot	P56501
GeneID	22229;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ucp3

Note: This product is for in vitro research use only