

Triosephosphate isomerase

Catalog No: #AP88841

Package Size: #AP88841-1 50ug #AP88841-2 100ug #AP88841-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Triosephosphate isomerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TIM,Triose-phosphate isomerase
Accession No.	Uniprot ID:Q5E956
Uniprot	Q5E956
GeneID	281543;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TPI1

Note: This product is for in vitro research use only