

Transcriptional coactivator YAP1

Catalog No: #AP88843



Package Size: #AP88843-1 50ug #AP88843-2 100ug #AP88843-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcriptional coactivator YAP1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Yes-associated protein 1,Protein yorkie homolog,Yes-associated protein YAP65 homolog,Yap,Yap65
Accession No.	Uniprot ID:Q2EJA0
Uniprot	Q2EJA0
GeneID	363014;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Yap1

Note: This product is for in vitro research use only