

## mRNA decay activator protein ZFP36

Catalog No: #AP88847



Package Size: #AP88847-1 50ug #AP88847-2 100ug #AP88847-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	mRNA decay activator protein ZFP36
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Growth factor-inducible nuclear protein NUP475,TPA-induced sequence 11,Tristetraprolin,Zinc finger protein 36,Zfp-36,Nup475,Tis11,Tis11a,Ttp
Accession No.	Uniprot ID:P22893
Uniprot	P22893
GeneID	22695;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Zfp36

Note: This product is for in vitro research use only