

## Sodium/glucose cotransporter 2

Catalog No: #AP88850



Package Size: #AP88850-1 50ug #AP88850-2 100ug #AP88850-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Sodium/glucose cotransporter 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Na(+)/glucose cotransporter 2,Low affinity sodium-glucose cotransporter,Solute carrier family 5 member 2,Sgt2
Accession No.	Uniprot ID:Q92317
Uniprot	Q92317
GeneID	246787;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Slc5a2

Note: This product is for in vitro research use only