

Tripartite motif-containing protein 3

Catalog No: #AP88852



Package Size: #AP88852-1 50ug #AP88852-2 100ug #AP88852-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tripartite motif-containing protein 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	RING finger protein 22,RING finger protein HAC1,Hac1,Rnf22
Accession No.	Uniprot ID:Q9R1R2
Uniprot	Q9R1R2
GeneID	55992;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Trim3

Note: This product is for in vitro research use only