

Trem-like transcript 1 protein

Catalog No: #AP88853



Package Size: #AP88853-1 50ug #AP88853-2 100ug #AP88853-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Trem-like transcript 1 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TLT-1, Triggering receptor expressed on myeloid cells-like protein 1, Tlt1
Accession No.	Uniprot ID:Q8K558
Uniprot	Q8K558
GeneID	71326;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Trem1

Note: This product is for in vitro research use only