

Tandem C2 domains nuclear protein

Catalog No: #AP88855



Package Size: #AP88855-1 50ug #AP88855-2 100ug #AP88855-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tandem C2 domains nuclear protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Membrane targeting tandem C2 domain-containing protein 1, Tandem C2 protein in nucleus, Tac2-N, C14orf47, MTAC2D1
Accession No.	Uniprot ID: Q8N9U0
Uniprot	Q8N9U0
GeneID	123036;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TC2N

Note: This product is for in vitro research use only