

Tumor necrosis factor receptor superfamily member 5

Catalog No: #AP88859



Package Size: #AP88859-1 50ug #AP88859-2 100ug #AP88859-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	B-cell surface antigen CD40,CD40L receptor,TNFRSF5,CD40
Accession No.	Uniprot ID:Q7YRL5
Uniprot	Q7YRL5
GeneID	403469;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD40

Note: This product is for in vitro research use only