

Tyrosine-protein kinase receptor UFO

Catalog No: #AP88860



Package Size: #AP88860-1 50ug #AP88860-2 100ug #AP88860-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tyrosine-protein kinase receptor UFO
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Adhesion-related kinase,Ark,Ufo
Accession No.	Uniprot ID:Q00993
Uniprot	Q00993
GeneID	26362;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Axl

Note: This product is for in vitro research use only