

Translational activator of cytochrome c oxidase 1

Catalog No: #AP88862



Package Size: #AP88862-1 50ug #AP88862-2 100ug #AP88862-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Translational activator of cytochrome c oxidase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Coiled-coil domain-containing protein 44, Translational activator of mitochondrially-encoded cytochrome c oxidase 1, Ccdc44
Accession No.	Uniprot ID:Q8K0Z7
Uniprot	Q8K0Z7
GeneID	70207;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Taco1

Note: This product is for in vitro research use only