

## Trimethyllysine dioxygenase, mitochondrial

Catalog No: #AP88864



Package Size: #AP88864-1 50ug #AP88864-2 100ug #AP88864-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Trimethyllysine dioxygenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Epsilon-trimethyllysine 2-oxoglutarate dioxygenase, TML hydroxylase, TML-alpha-ketoglutarate dioxygenase, TML dioxygenase
Accession No.	Uniprot ID:Q0VC74
Uniprot	Q0VC74
GeneID	535630;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: TMLHE

Note: This product is for in vitro research use only