

Tumor necrosis factor receptor superfamily member 4

Catalog No: #AP88867



Package Size: #AP88867-1 50ug #AP88867-2 100ug #AP88867-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	MRC OX40,OX40 antigen,OX40L receptor,CD134,Ox40,Txgp11
Accession No.	Uniprot ID:P15725
Uniprot	P15725
GeneID	25572;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnfrsf4

Note: This product is for in vitro research use only