

Vascular endothelial growth factor receptor 1

Catalog No: #AP88872



Package Size: #AP88872-1 50ug #AP88872-2 100ug #AP88872-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vascular endothelial growth factor receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	VEGFR-1,Embryonic receptor kinase 2,Fms-like tyrosine kinase 1,FLT-1,Tyrosine-protein kinase receptor FLT,Emrk2,Flt,Vegfr1
Accession No.	Uniprot ID:P35969
Uniprot	P35969
GeneID	14254;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Flt1

Note: This product is for in vitro research use only