

## TGF-beta receptor type-1

Catalog No: #AP88874



Package Size: #AP88874-1 50ug #AP88874-2 100ug #AP88874-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TGF-beta receptor type-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TGFR-1,Activin A receptor type II-like protein kinase of 53kD,Activin receptor-like kinase 5,ALK-5,Serine/threonine-protein kinase receptor R4,SKR4,TGF-beta type I receptor,Transforming growth factor-beta receptor type I,TGF-beta receptor type I,ALK5,SKR4
Accession No.	Uniprot ID:P36897
Uniprot	P36897
GeneID	7046;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TGFBR1

Note: This product is for in vitro research use only