

Tsukushin

Catalog No: #AP88878



Package Size: #AP88878-1 50ug #AP88878-2 100ug #AP88878-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tsukushin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tsukushi,Leucine-rich repeat-containing protein 54,Lrrc54,Tsk
Accession No.	Uniprot ID:Q8CBR6
Uniprot	Q8CBR6
GeneID	244152;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tsku

Note: This product is for in vitro research use only