

Transthyretin

Catalog No: #AP88879



Package Size: #AP88879-1 50ug #AP88879-2 100ug #AP88879-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transthyretin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Prealbumin,TBPA,Tt
Accession No.	Uniprot ID:P02767
Uniprot	P02767
GeneID	24856;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ttr

Note: This product is for in vitro research use only