

Toll-interacting protein

Catalog No: #AP88884

Package Size: #AP88884-1 50ug #AP88884-2 100ug #AP88884-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Toll-interacting protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	RCJMB04_14a3
Accession No.	Uniprot ID:Q5ZK05
Uniprot	Q5ZK05
GeneID	423099;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TOLLIP

Note: This product is for in vitro research use only