

Tryptophan 5-hydroxylase 1

Catalog No: #AP88885



Package Size: #AP88885-1 50ug #AP88885-2 100ug #AP88885-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tryptophan 5-hydroxylase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Tryptophan 5-monoxygenase 1,TPH
Accession No.	Uniprot ID:P70080
Uniprot	P70080
GeneID	395799;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TPH1

Note: This product is for in vitro research use only