

Toll/interleukin-1 receptor domain-containing adapter protein

Catalog No: #AP88888

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP88888-1 50ug #AP88888-2 100ug #AP88888-3 1mg

Description

Product Name	Toll/interleukin-1 receptor domain-containing adapter protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TIR domain-containing adapter protein,Adaptor protein Wyatt,MyD88 adapter-like protein,Mal
Accession No.	Uniprot ID:Q99JY1
Uniprot	Q99JY1
GeneID	117149;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tirap

Note: This product is for in vitro research use only