

Troponin C, slow skeletal and cardiac muscles

Catalog No: #AP88889



Package Size: #AP88889-1 50ug #AP88889-2 100ug #AP88889-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Troponin C, slow skeletal and cardiac muscles
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TN-C,Tncc
Accession No.	Uniprot ID:P19123
Uniprot	P19123
GeneID	21924;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnnc1

Note: This product is for in vitro research use only