

Troponin T, fast skeletal muscle

Catalog No: #AP88893



Package Size: #AP88893-1 50ug #AP88893-2 100ug #AP88893-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Troponin T, fast skeletal muscle
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TnTf, Fast skeletal muscle troponin T, fTnT
Accession No.	Uniprot ID:Q9QZ47
Uniprot	Q9QZ47
GeneID	21957;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tnnt3

Note: This product is for in vitro research use only