

Tumor necrosis factor receptor superfamily member 16

Catalog No: #AP88894

Package Size: #AP88894-1 50ug #AP88894-2 100ug #AP88894-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 16
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Gp80-LNGFR,Low affinity neurotrophin receptor p75NTR,Low-affinity nerve growth factor receptor,NGF receptor,p75 ICD,CD271,TNFRSF16
Accession No.	Uniprot ID:P08138
Uniprot	P08138
GeneID	4804;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NGFR

Note: This product is for in vitro research use only