

Fermitin family homolog 3

Catalog No: #AP88898



Package Size: #AP88898-1 50ug #AP88898-2 100ug #AP88898-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fermitin family homolog 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Kindlin-3,MIG2-like protein,Unc-112-related protein 2,KIND3,MIG2B,URP2
Accession No.	Uniprot ID:Q86UX7
Uniprot	Q86UX7
GeneID	83706;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:FERMT3

Note: This product is for in vitro research use only