

Protein unc-13 homolog C

Catalog No: #AP88900



Package Size: #AP88900-1 50ug #AP88900-2 100ug #AP88900-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein unc-13 homolog C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Munc13-3,Unc13h3
Accession No.	Uniprot ID:Q62770
Uniprot	Q62770
GeneID	286931;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Unc13c

Note: This product is for in vitro research use only