

RISC-loading complex subunit TARBP2

Catalog No: #AP88902



Package Size: #AP88902-1 50ug #AP88902-2 100ug #AP88902-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	RISC-loading complex subunit TARBP2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TAR RNA-binding protein 2, Trans-activation-responsive RNA-binding protein, TRBP
Accession No.	Uniprot ID:Q15633
Uniprot	Q15633
GeneID	6895;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TARBP2

Note: This product is for in vitro research use only