

# Tubulin alpha-1A chain

Catalog No: #AP88906



Package Size: #AP88906-1 50ug #AP88906-2 100ug #AP88906-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tubulin alpha-1A chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alpha-tubulin 1, Tubulin alpha-1 chain, Tuba1
Accession No.	Uniprot ID:P68370
Uniprot	P68370
GeneID	64158;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Tuba1a

Note: This product is for in vitro research use only