

Vasopressin V1b receptor

Catalog No: #AP88908



Package Size: #AP88908-1 50ug #AP88908-2 100ug #AP88908-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vasopressin V1b receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	V1bR,AVPR V1b,AVPR V3,Antidiuretic hormone receptor 1b,Vasopressin V3 receptor
Accession No.	Uniprot ID:P48974
Uniprot	P48974
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Avpr1b

Note: This product is for in vitro research use only