

GATA-type zinc finger protein 1

Catalog No: #AP88913

Package Size: #AP88913-1 50ug #AP88913-2 100ug #AP88913-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | GATA-type zinc finger protein 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Ms |
| Other Names | GATA-like protein 1, GLP-1, Glp1 |
| Accession No. | Uniprot ID:Q1WG82 |
| Uniprot | Q1WG82 |
| GeneID | 100009600; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Zglp1

Note: This product is for in vitro research use only