

T-box transcription factor TBX3

Catalog No: #AP88914



Package Size: #AP88914-1 50ug #AP88914-2 100ug #AP88914-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-box transcription factor TBX3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	T-box protein 3
Accession No.	Uniprot ID:Q7TST9
Uniprot	Q7TST9
GeneID	353305;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tbx3

Note: This product is for in vitro research use only