

# TGF-beta-activated kinase 1 and MAP3K7-binding protein 1

Catalog No: #AP88915

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP88915-1 50ug #AP88915-2 100ug #AP88915-3 1mg

## Description

Product Name	TGF-beta-activated kinase 1 and MAP3K7-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mitogen-activated protein kinase kinase kinase 7-interacting protein 1,TGF-beta-activated kinase 1-binding protein 1,TAK1-binding protein 1,Map3k7ip1
Accession No.	Uniprot ID:Q8CF89
Uniprot	Q8CF89
GeneID	66513;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tab1

Note: This product is for in vitro research use only