

Tapasin-related protein

Catalog No: #AP88917



Package Size: #AP88917-1 50ug #AP88917-2 100ug #AP88917-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tapasin-related protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TAPASIN-R, TAP-binding protein-like, TAP-binding protein-related protein, TAPBP-R, Tapasin-like
Accession No.	Uniprot ID:Q9BX59
Uniprot	Q9BX59
GeneID	55080;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TAPBPL

Note: This product is for in vitro research use only