

TYRO protein tyrosine kinase-binding protein

Catalog No: #AP88919



Package Size: #AP88919-1 50ug #AP88919-2 100ug #AP88919-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	TYRO protein tyrosine kinase-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	DNAX-activation protein 12,Dap12
Accession No.	Uniprot ID:Q6X9T7
Uniprot	Q6X9T7
GeneID	361537;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tyrobp

Note: This product is for in vitro research use only