

L-dopachrome tautomerase

Catalog No: #AP88920



Package Size: #AP88920-1 50ug #AP88920-2 100ug #AP88920-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-dopachrome tautomerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	DCT,L-dopachrome Delta-isomerase,Tyrosinase-related protein 2,TRP-2,TYRP2
Accession No.	Uniprot ID:O93505
Uniprot	O93505
GeneID	395775;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DCT

Note: This product is for in vitro research use only