

Thyrotropin subunit beta

Catalog No: #AP88922



Package Size: #AP88922-1 50ug #AP88922-2 100ug #AP88922-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thyrotropin subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Thyroid-stimulating hormone subunit beta, TSH-B, Thyrotropin beta chain
Accession No.	Uniprot ID:O57340
Uniprot	O57340
GeneID	395937;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TSHB

Note: This product is for in vitro research use only