

Thioredoxin reductase 1, cytoplasmic

Catalog No: #AP88925



Package Size: #AP88925-1 50ug #AP88925-2 100ug #AP88925-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thioredoxin reductase 1, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TR, Gene associated with retinoic and interferon-induced mortality 12 protein, GRIM-12, KM-102-derived reductase-like factor, Thioredoxin reductase TR1, GRIM12, KDRF
Accession No.	Uniprot ID:Q16881
Uniprot	Q16881
GeneID	7296;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TXNRD1

Note: This product is for in vitro research use only