

Anionic trypsin-1

Catalog No: #AP88926



Package Size: #AP88926-1 50ug #AP88926-2 100ug #AP88926-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Anionic trypsin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Anionic trypsin I,Pretrypsinogen I,Serine protease 1,Try1
Accession No.	Uniprot ID:P00762
Uniprot	P00762
GeneID	24691;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prss1

Note: This product is for in vitro research use only