

Urotensin-2B

Catalog No: #AP88929



Package Size: #AP88929-1 50ug #AP88929-2 100ug #AP88929-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Urotensin-2B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Urotensin II-related peptide,Urotensin IIB,U-IIB,Urotensin-2 domain-containing protein,Urp,Uts2d
Accession No.	Uniprot ID:Q76512
Uniprot	Q76512
GeneID	378939;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Uts2b

Note: This product is for in vitro research use only