

Serp H1

Catalog No: #AP88930



Package Size: #AP88930-1 50ug #AP88930-2 100ug #AP88930-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serp H1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	47 kDa heat shock protein, Collagen-binding protein, Colligin, Serine protease inhibitor J6, Cbp1, Hsp47
Accession No.	Uniprot ID:P19324
Uniprot	P19324
GeneID	12406;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Serpinh1

Note: This product is for in vitro research use only