

Sestrin-2

Catalog No: #AP88931



Package Size: #AP88931-1 50ug #AP88931-2 100ug #AP88931-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sestrin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Hypoxia-induced gene,HI95,SEST2
Accession No.	Uniprot ID:P58004
Uniprot	P58004
GeneID	83667;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SESN2

Note: This product is for in vitro research use only