

Urea transporter 1

Catalog No: #AP88932



Package Size: #AP88932-1 50ug #AP88932-2 100ug #AP88932-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Urea transporter 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Solute carrier family 14 member 1,Urea transporter B,UT-B,Urea transporter, erythrocyte,UT-B
Accession No.	Uniprot ID:Q5QF96
Uniprot	Q5QF96
GeneID	493988;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SLC14A1

Note: This product is for in vitro research use only