

Thiopurine S-methyltransferase

Catalog No: #AP88933



Package Size: #AP88933-1 50ug #AP88933-2 100ug #AP88933-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thiopurine S-methyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Thiopurine methyltransferase
Accession No.	Uniprot ID:Q8HX86
Uniprot	Q8HX86
GeneID	403536;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TPMT

Note: This product is for in vitro research use only