

Thrombopoietin

Catalog No: #AP88934



Package Size: #AP88934-1 50ug #AP88934-2 100ug #AP88934-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Thrombopoietin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C-mpl ligand,ML,Megakaryocyte colony-stimulating factor,Megakaryocyte growth and development factor,MGDF,Myeloproliferative leukemia virus oncogene ligand
Accession No.	Uniprot ID:P40226
Uniprot	P40226
GeneID	21832;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Thpo

Note: This product is for in vitro research use only