

TNF receptor-associated factor 4

Catalog No: #AP88938



Package Size: #AP88938-1 50ug #AP88938-2 100ug #AP88938-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TNF receptor-associated factor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cysteine-rich motif associated to RING and Traf domains protein 1, Cart1
Accession No.	Uniprot ID:Q61382
Uniprot	Q61382
GeneID	22032;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Traf4

Note: This product is for in vitro research use only