

Tumor necrosis factor ligand superfamily member 18

Catalog No: #AP88939



Package Size: #AP88939-1 50ug #AP88939-2 100ug #AP88939-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 18
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Activation-inducible TNF-related ligand,AITRL,Glucocorticoid-induced TNF-related ligand,hGITRL,AITRL,GITRL,TL6,UNQ149/PRO175
Accession No.	Uniprot ID:Q9UNG2
Uniprot	Q9UNG2
GeneID	8995;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNFSF18

Note: This product is for in vitro research use only