

Human immunodeficiency virus type I enhancer-binding protein 2 homolog

Catalog No: #AP88942

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP88942-1 50ug #AP88942-2 100ug #AP88942-3 1mg

Description

| | |
|--------------------|---|
| Product Name | Human immunodeficiency virus type I enhancer-binding protein 2 homolog |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | Angiotensinogen gene-inducible enhancer-binding protein 1,DNA-binding protein AGIE-BP1,Myc intron-binding protein 1,MIBP-1,Agie-bp1,Mibp1 |
| Accession No. | Uniprot ID:Q00900 |
| Uniprot | Q00900 |
| GeneID | 29721; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hivep2

Note: This product is for in vitro research use only