

WW domain-containing oxidoreductase

Catalog No: #AP88943



Package Size: #AP88943-1 50ug #AP88943-2 100ug #AP88943-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WW domain-containing oxidoreductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Fragile site FRA16D oxidoreductase, Short chain dehydrogenase/reductase family 41C member 1, FOR, SDR41C1, WOX1
Accession No.	Uniprot ID: Q9NZC7
Uniprot	Q9NZC7
GeneID	51741;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: WWOX

Note: This product is for in vitro research use only