

WW domain-containing oxidoreductase

Catalog No: #AP88943

Package Size: #AP88943-1 50ug #AP88943-2 100ug #AP88943-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	WW domain-containing oxidoreductase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Fragile site FRA16D oxidoreductase,Short chain dehydrogenase/reductase family 41C member 1,FOR,SDR41C1,WOX1
Accession No.	Uniprot ID:Q9NZC7
Uniprot	Q9NZC7
GeneID	51741;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:WVOX

Note: This product is for in vitro research use only